Transcutaneous Pacing Procedure

- See 14. Adult Bradycardia (< 50/min)/Tachycardia (> 150/min) (with Pulse) for indications. Inotropes may be used if transcutaneous pacing is NOT available. See 14. Adult Bradycardia (< 50/min)/Tachycardia (> 150/min) (with Pulse)
- 2. Place the pacing pads on the chest of the patient as per package instructions
- 3. Connect the pads cable to the pacing machine if not already connected
- **4. Turn the pacer ON**. Observe for **markers** (•) indicating the R-wave on the screen. Some machines require that you **START pacing** after turning the pacer on. Observe for **pacing spikes** (|) on the baseline.
- **5.** Set the **Rate** to approximately 60-70 bpm.
- 6. Set current milliamperes (mA) output as follows: Increase milliamperes (mA) from minimum setting until every pacer spike is immediately followed by a wide QRS and a broad T wave This is termed as Electrical Capture.
- 7. Confirm by checking the patients femoral pulse to see if the pulse rate matches the rate set above i.e. 60-70bpm. This is termed as **Mechanical Capture.**
- 8. Recheck the patient's vital signs and confirm the patient's signs of shock are resolving i.e. increase in blood pressure, improved mentation, etc. This is termed as Physiological Capture.
- 9. If all the above is achieved, increase the current milliamperes by 10% for safety margin
- **10.** Provide adequate sedation and analgesia if the patient experiences any discomfort
- Transfer care to a Cardiologist without delay. DO NOT STOP PACING unless instructed to by a Cardiologist.

Trouble Shooting

- Pacing Spikes not seen on the base line Confirm that you have pressed the START button
- **No Electrical Capture** Confirm that the pads are firmly pressed on the patient's chest. Continue increasing the milliamperes. There is no set minimum or maximum.
- **No Mechanical Capture** Increase the milliamperes by increments of 5-10mA and recheck the pulse
- No Physiological Capture Consider hypovolaemia as the cause of shock and give a small fluid bolus (250-500mls) and recheck the patient. If not, increase the set rate to 80bpm, confirm electrical capture and mechanical capture and recheck the patient
- In all cases, consult a Cardiologist.